

# Viral RNA/DNA Purification Kit

(Spin Columns)

Qualitative Assay for  
Manual Extraction Systems

## Instructions For Use



Instructions For Use  
Version 2. GE-009.03.23



GE-009

GE\_009/50 – Viral RNA/DNA Purification Kit – 50 rxn

GE\_009/250 – Viral RNA/DNA Purification Kit – 250 rxn

GE\_009/500 – Viral RNA/DNA Purification Kit – 500 rxn



OxGEN Solutions, 14<sup>th</sup> km, Natakhtari,  
Mtskheta, Municipality 3308, Georgia



MDSS GmbH 41 Schiffgraben,  
30175, Hannover, Germany

## Description of the Kit Components, Transportation and Storage

Table 1. Content of the Kit

Viral RNA/DNA Purification Kit		50 rxn	250 rxn	500 rxn	Transportation and Storage
Solution A	Lysis/Binding Buffer	15 ml	71 ml	71 ml x 2	Room temperature
Solution W1 (conc.)	Wash Buffer 1	8 ml	21 ml x 2	28 ml x 3	Room temperature
Solution W2 (conc.)	Wash Buffer 2	5 ml	14 ml x 2	18 ml x 3	Room temperature
Solution E	Elution Buffer	11 ml	51 ml	105 ml	Room temperature
Carrier RNA	Synthetic Oligo Mix <i>Lyophilized Powder</i>	25 µg	125 µg	125 µg x 2	Room temperature
G-Spin/Columns	Silica Spin Columns, with Collections	50	250	500	Room temperature
Collection Tubes	Collection Tubes (2 ml)	100	500	1000	Room temperature

## Reagents Preparation

### Solution W1

Wash Buffer 1 comes as a concentrate. Prior to initial use, combine the recommended quantity of ethanol, which must be at least 95% pure, as specified on the bottle and in Table 2. If the labels of Solution W1 indicate that ethanol has already been added by the manufacturer  omit this step.

Table 2. Preparation of Solution Wash 1

No. Reactions	Solution W1	Ethanol ≥95%	Final Volume
50	8 ml	23 ml	31 ml
250	21 ml x 2	59 ml (In each bottle)	80 ml x 2 bottles
500	28 ml x 3	77 ml (In each bottle)	105 ml x 3 bottles

## Solution W2

Wash Buffer 2 comes as a concentrate. Prior to initial use, combine the recommended quantity of ethanol, which must be at least 95% pure, as specified on the bottle and in Table 3. If the labels of Solution W2 indicate that ethanol has already been added by the manufacturer  omit this step.

Table 3. Preparation of Solution Wash 2

No. Reactions	Solution W2	Ethanol $\geq 95\%$	Final Volume
50	5 ml	26 ml	31 ml
250	14 ml x 2	66 ml (In each bottle)	80 ml x 2 bottles
500	18 ml x 3	87 ml (In each bottle)	105 ml x 3 bottles

## Carrier RNA

Before using it for the first time, add the appropriate amount of Solution E (Elution Buffer) as indicated in Table 4. Aliquoted Carrier RNA should be stored at  $-20^{\circ}\text{C}$ . Do not freeze-thaw more than 5 times.

Table 4. Preparation of Carrier RNA

No. Reactions	Carrier RNA	Solution E	Final Volume
50	25 $\mu\text{g}$	0.25 ml	0.25 ml x 1 vial
250	125 $\mu\text{g}$	1.25 ml	1.25 ml x 1 vials
500	125 $\mu\text{g}$ x 2	1.25 ml (In each vial)	1.25 ml x 2 vials

## Recommended Sample Pretreatment

This kit facilitates the extraction of RNA and DNA from viruses.

**Urine** - Please centrifuge 2 ml urine samples at 2,250 g for 15 minutes before processing. After centrifugation, resuspend the resulting pellet in 200 µl of 1x PBS buffer.

**Swabs (Nasopharyngeal, Oropharyngeal, Urogenital)** - Swab samples may contain mucus, blood, or other materials that can negatively affect purification. To prepare for processing, thoroughly mix the swab in 1x PBS or 0.9 % saline buffer. Vortex the mixture vigorously at the highest setting for 15 seconds, then transfer the clear sample into a microwell plate.

**Whole Blood** - Before processing, ensure that entire blood samples are fully liquefied to prevent clot carryover, which could interfere with nucleic acid purification.

**Plasma/Serum** - When handling plasma and serum samples, carefully transfer the clear, yellowish fluid, not disturbing the underlying buffy coat or red blood cell layer.

**Semen/Cerebrospinal Fluid** - Vortex the samples vigorously for no more than 15 seconds, and then transfer the clarified sample into a microwell plate or microcentrifuge tube.

## The List of Materials to be Supplied by the User

Table 5. Equipment and Reagents to be Supplied by the User

Equipment	Consumables
Thermoblock or thermomixer	Ethanol ≥95%
Centrifuge and minifuge	RNase-free 1.5 ml microcentrifuge tubes
Vortex	0.5 - 10 µl pipette tips with filter
Pipette 0.5 - 10 µl	10 – 100 µl pipette tips with filter
Pipette 20 - 200 µl	20 – 200 µl pipette tips with filter
Pipette 100 - 1000 µl	100 - 1000 µl pipette tips with filter

## Instructions for Manual Purifications

**IMPORTANT:** The sample should be stored according to “Collecting and Handling of Clinical Specimens for PCR Testing”. Before analysis, ensure the samples do not contain inhibiting mucus and sediments.

### Viral RNA/DNA Purification Kit (Spin Columns) Protocol

**Note:** Before starting the procedure, prepare the solutions according to the solution preparation guide (Table 2-4). Solution A may form precipitates upon storage. Warm it up to 60°C until the residues have fully dissolved. The mixture of Solution A and Carrier RNA can be prepared in advance and is stable at 2-8°C for up to 72 h (for one reaction, add 5 µl of Carrier RNA to 280 µl of Solution A). **Preheat Solution E at 56°C before starting the procedure.**

1. Transfer 200 µl of sample into a 1.5 ml tube, add 280 µl of Solution A and 5 µl of Carrier RNA. Mix by pulse-vortexing for 15 sec;  
**Note:** If the Solution A and Carrier RNA mix are prepared in advance, transfer 285 µl of the mix into the sample tube.
2. Incubate the sample for 15 min at 56°C in a thermomixer at 1400 rpm. Alternatively, incubate in a thermoblock and vortex periodically at 5 min intervals and spin down;
3. Cool down the sample at Room Temperature (RT). Add 200 µl of 96–100% cold ethanol, invert gently 30 times and spin down.
4. Transfer 680 µl of lysate into a G-spin/column and centrifuge at 8 000 rpm for 2 min. Discard the collection tube;
5. Wash the G-spin/column with 600 µl of Solution W1 and centrifuge at 13 000 rpm for 2 min. Discard the collection tube;
6. Wash the G-spin/column with 600 µl of Solution W2 centrifuge at 13 000 rpm for 2 min. Discard the flow-through;
7. Remove residual wash buffer by centrifuging at 13 000 rpm for 2 min. Discard the collection tube;
8. Transfer the G-spin/column into a new 1.5 ml microcentrifuge tube;
9. Add 50-200 µl of preheated (56°C) Solution E to the G-spin/column, ensuring the membrane's entire surface is hydrated. Avoid touching the membrane with the pipette tip;
10. Incubate the samples for 3 min at RT.
11. Elute the RNA/DNA by centrifuging at 13 000 rpm for 1 min.

## Disposal

Dispose of used kit reagents, human clinical samples, and sealed amplification plates as laboratory clinical waste according to local, state, and federal regulations.

## Version History

Instruction for Use Version GE-009.03.23 EN V2, March 10, 2023.

## Quality Control System

Quality management system TÜV SÜD-ISO 9001:2015. Each Viral RNA/DNA Purification Kit batch is tested against predetermined quality specifications to ensure consistent product quality.

## Technical support

For technical support, please contact our dedicated Technical Support Team at:  
TEL: +995 599 374 374, Email: [support@oxgensolutions.com](mailto:support@oxgensolutions.com)

## Explanation of Symbols

**LOT** Batch code

 Use by

**REF** Catalogue number

 Store at

**QTY** Quantity

 Manufactured by